

SPECIFICATIONS

tCOM+ with V-Sign™ Sensor 2

The tCOM+ Transcutaneous Monitor features an intuitive interface that displays parameters as values, trends, and deltas with the ability to set preferred alarm limits with onboard tutorials that support smooth caregiver workflows.

The digital V-Sign™ sensor enables continuous and noninvasive monitoring of tcPCO₂, SpO₂, PR (pulse rate), and RHP (relative heating power).

Suitable for Adults, Pediatrics & Neonates. SpO₂ and PR available for patients aged term birth plus 12 months.



Overall System Performance

Transcutaneous Carbon Dioxide Partial Pressure (tcPCO₂)¹

Measurement range: 0–200 mmHg (0–26.7 kPa)

Resolution: 0.1 mmHg (0.01 kPa) below 100 mmHg (10 kPa)/1 mmHg (0.1 kPa) above 100 mmHg (10 kPa)

Drift²: Typically < 0.5%/hour

Response time (T90)²: Typically < 75 sec

Linearity²: Typically < 1 mmHg (0.13 kPa)

Interferences by anesthetic gases: Negligible

Stabilization/artifact detection: After sensor application or occurrence of a tcPCO₂ artifact, tcPCO₂ is displayed in grey until it (re)stabilizes.

Approved sites for tcPCO₂ monitoring: Monitoring in neonates include the thorax (under the clavicle, on rib cage), abdomen, back, low forehead, and the inner or anterior aspect of the thigh. Approved sites for PCO₂ monitoring in adult and pediatric patients include the earlobe, low forehead, cheek, upper arm, scapula (shoulder blade), thorax (under the clavicle), mastoid process, and lower arm.

Pulse Rate (PR)

Available for patients older than term birth + 12 months.

Measurement range: 30–250 bpm

Resolution: 1 bpm

Accuracy: ± 3 bpm

Pulsation Index (PI)

Available for patients older than term birth + 12 months.

Measurement range: 0.1–10.0%

Resolution: 0.1%

Oxygen Saturation (SpO₂)

Available for patients older than term birth + 12 months.

Measurement range: 1–100%

Resolution: 1%

Accuracy: ±2% (A_{rms} over 70% to 100%)³

Averaging mode: 2, 3, 4, 6, 8, 12, 16, and 32 sec

Approved sites for SpO₂/ PR monitoring: Earlobe, low forehead, cheek, upper arm, and scapula (shoulder blade)

Sensor Heating Power (HP)

Measurement range:

Absolute Heating Power (AHP): 0–999 mW

Relative Heating Power (RHP): -999–999 mW

Resolution: 1 mW

Sensor Temperature

Measurement range: 0.0–70.0°C

Resolution: 0.1°C

Accuracy: ± 0.2°C (over 37.0 to 45.0°C)

¹ An algorithm developed by J.W. Severinghaus is used to calculate tcPCO₂ from the measured cutaneous PCO₂. This algorithm accounts for temperature and metabolic correction factors. The tcPCO₂ values displayed by the tCOM+ are corrected / normalized to 37°C and provide an estimate of arterial PCO₂ (PaCO₂) at 37°C. Correction factors can be customized by institution. Additionally, and subject to institution's permission, In-vivo Correction (IIC) of tcPCO₂ values is possible at the bedside.

² Respective specifications based on in vitro tests performed as per IEC 60601-2-23:2011 at a sensor temperature of 42°C.

³ SpO₂ accuracy specification is based on controlled hypoxia studies on healthy, adult volunteers over the specified saturation range by applying a V-Sign™ Sensor 2 to each of the specified measurement sites.

V-Sign™ Sensor 2 (Digital, combined tcPCO₂/ SpO₂/ PR sensor)



General Characteristics

Suitable for adult, pediatric, and neonatal patients
Reusable, waterproof

Measurement Principle

Electrochemical Severinghaus-type tcPCO₂ sensor combined with reflectance 2-wavelength pulse oximetry.

Digital Microtechnology

Highly integrated opto-electronic sensor head comprising micro pH-electrode, optical oximetry unit, temperature sensors, heating unit, all combined in a fully digital design. High definition digitizer and pre-processing in the sensor head provides robust and low noise signals that are digitally transmitted to the tCOM+.

Sensor Memory

Sensor-specific data are stored in the sensor's memory after manufacturing (serial number, factory PCO₂ sensitivity / calibration, etc.) and during operation (sensor calibration, membrane change, etc.).

Sensor Membrane Change

The default membrane change interval is 28 days but may be customized depending on the specific requirements of various clinical settings. The membrane changer tool simplifies the process for a highly reproducible membrane change.

Sensor Calibration

Calibration duration: Typically 3 minutes (ex. factory)

Calibration interval: Up to 12 hours. Once the calibration interval has elapsed, sensor calibration is **recommended**, monitoring is possible for another 4 to 6 hours with tcPCO₂ marked as 'questionable'. Thereafter, sensor calibration is **mandatory** and tcPCO₂ is marked as 'invalid' (tcPCO₂ value replaced by '---').

Smart Cal-Mem

Supports the disconnection of the sensor for up to 30 minutes without losing the calibration status. Furthermore, the sensor can be removed from the Docking Station for up to 10 minutes without initiating a calibration upon reinsertion of the sensor into the Docking Station. Overall, Smart Cal-Mem significantly reduces the number of required calibrations and the calibration gas consumption.

Sensor Internal Temperature Control

To ensure safe operation, the sensor reliably supervises the sensor temperature with two independent circuits. Additionally, the tCOM+ software redundantly controls the temperature of the connected sensor for an extra level of protection.

Sensor Dimensions / Sensor Cable

Diameter x height of sensor head: 14 mm x 9 mm (0.55" x 0.35")

Weight of sensor head: < 2.9 g (0.1 oz)

Length of sensor cable: Approximately 80 cm (31"), plus Sensor Adapter Cable available in lengths of Regular (150 cm (59")), Long (250 cm (98")), or Extra-Long (750 cm (295")).

Sensor cable: Highly flexible, shielded cable with coating withstanding cleaning agents and irradiation commonly used in busy hospital environments.

Transport / Storage of Sensor

Transport temperature: 0–50°C (32–122°F)

Long term storage temperature: 15–26°C (59–78°F)

Transport and store sensor with membrane applied; protect sensor from light and radiation.

tCOM+ Monitor

Physical Characteristics

Weight: 2.5 kg (5.5 lbs)

Size: 15.3 cm (height) x 27.8 cm (width) x 16.2 cm (depth) (6.02" x 10.95" x 6.38")

Transport: Foldable handle located on top of the monitor for monitor transport on the unit.

Mounting: Compatible with 75x75 VESA mounts for roll/infusion stands, wall / railing mounts, transport incubators, etc.

Cable storage: Included cable holder may be attached to the back right or left of the monitor to stow the cable during transport or storage.

Tilting: Option to add feet to the VESA holes to adjust angle for improved table-top viewing.

Sensor Calibration

Built-in sensor calibration chamber for 1-point calibration. Automatic calibration ensures that system is 'Ready for use' if sensor is stored in calibration chamber. Comprehensive controls ensure reliable calibrations.

Sensor Temperature

Selectable range: Configurable by institution between 37.0 and 44°C (Neonate) and 37.0 and 44.5°C (Adult/Pediatric). Safety controls of the monitor may restrict the selectable range depending on the type of the connected sensor, the selected patient mode, or the enabled parameters.

Default sensor temperature: 42.0°C in Adult Mode and 41.0°C in Neonatal Mode.

Site Time

Selectable 'Site Time': Configurable by the institution between 0.5 and 12.0 hours maximum (in steps of 0.5 hours). Depending on the selected patient mode and with increasing sensor temperature, safety controls of the monitor may enforce a shorter site time.

Default 'Site Time': 8.0 hours in Adult Mode at 42.0°C or in Neonatal Mode at 41.0°C (or as restricted by institution and / or safety controls of the monitor).

Site Timer

Timer indicating remaining 'Site Time' during monitoring. An alarm is triggered once 'Site Time' has elapsed.

Site Protection

Safety feature which prevents excessively long exposure of the skin to temperatures exceeding 41°C in adult mode or 40°C in neonatal mode. When Site Protection is active and 'Site Time' has elapsed, the sensor temperature will be reduced to 39°C for Neonatal and Adult Mode (41°C if SpO₂ enabled in Adult Mode).

Initial Heating

Temporarily increases sensor temperature after sensor application for faster perfusion and results (in Adult Mode only; +2°C max. 44.5°C; can only be switched on if enabled by institution).

Redundant Sensor Temperature Control on tCOM+

Should the sensor's temperature control fail, the monitor firmware also controls the temperature of the connected sensor. In the event of an error, the sensor will restart or switch off.

Alarm System

Alarm signals: Visual/ auditory alarm signals for high / low tcPCO₂, SpO₂, PR, and technical alarms.

Alarm inhibition: Auditory alarm signals can be PAUSED (1 or 2 minutes) or switched off permanently (if enabled by institution).

Alarm system status indicators: The on-screen Alarm Bar and LED Bar indicate the active alarm with the highest priority. The highest priority alarm is always displayed in the Alarm Bar. Auditory alarms can be paused for 1 or 2 minutes, or permanently switched off with the AUDIO OFF Reminder (unless disabled by the institution. Visual and auditory alarms may be turned off if option is enabled by the institution.

Display/ Indicators

LED Bar: Multiple colors to indicate high, medium, and low priority alarms

Display: 8.8" wide-angle, high color (16.7M) IPS display

Data update rate: 1 sec for tcPCO₂, SpO₂, PR, RHP

Data validity: Clear representation of data validity / quality for tcPCO₂, SpO₂, PR, PI, RHP (valid, questionable, unstable, invalid)

Measurement screens: Various preconfigured, user-selectable measurement screens displaying values and real time trends for tcPCO₂, SpO₂, PR, RHP; alarm limits for tcPCO₂, SpO₂, PR; Baseline, Baseline values and Delta-x values for tcPCO₂, SpO₂, RHP; values for pulsation index, AHP; Pleth Wave or blip bar reflecting relative pulse amplitude; visual alarm signals; status icons (e.g., remaining monitoring time) and status messages.

'Calibration' / 'Ready for use' screens: Display important system information about patient mode, profile, and calibration status and/or progress.

Tutorials: Tutorials provide step-by-step instructions for common workflows, including monitoring-specific tutorials like sensor application and sensor cleaning and maintenance-related tutorials like membrane change, sensor clean and soak, and sensitivity test. Certain troubleshooting workflows have interactive guidance available to help users resolve errors in real-time.

Favorites: Users may add different monitor menus to their Favorites tab to allow easy access to the most frequently used monitor menus.

Events: Users have the option to register events, i.e., medication administration, manipulation, or another user-created event. Several events can be stored in the tCOM+ for subsequent display on the measurement screen or in the Events log.

Languages: Catalan, Danish, Dutch, English, Finnish, French, German, Italian, Norwegian, Polish, Portuguese, Spanish, Swedish

Highly configurable: Patient Mode, Enabled Parameters, Severinghaus Correction Mode, Heating Power Mode, Spot Check Mode, Parameter Display Color, PCO₂/PO₂ Unit (mmHg or kPa), time ranges for online trends, sweep speed of Pleth wave, Sleep Mode, Brightness, Audio Settings, Menu Access, etc.

tCOM+ Profiles

Profiles help to ensure standardization of settings on all the monitors in a unit. Profiles can be configured in the password protected Advanced Settings menu to fit the specific needs of the clinical setting.

Special (safety-relevant) tCOM+ Parameters

Within a password-protected area of the monitor, the institution can configure all menu accessible parameters and special parameters to disable or restrict operator access to menu accessible parameters. For example, the maximal Sensor Temperature or Site Time can be adapted and restricted to safe settings for the patient population.

Patient Data Management

The tCOM+ automatically stores PCO₂, SpO₂, PR, RHP, and PI data in the internal memory for on-screen viewing of trends (up to 48 hours) and data export via the DATA/SERVICE USB port. The Data Recording Interval is 1 second; non-volatile memory provides storage of up to one year.

Data Communication

Wired

- Serial Data Port (RS-232): For communication with external data collection systems such as personal computers or multi-parameter bedside monitors.
- LAN port: for communication with external computer-based data collection systems.
- Analog output (0-1V): provides voltage outputs for tcPCO₂, SpO₂, PR, and the Pleth waveform (selectable ranges).
- DATA/SERVICE USB Port: for SW upgrades or data download via USB-C.
- Universal Connectivity Port: standard USB-C port for future connectivity options.

Wireless

- Wi-Fi: for connection to external networks and communication with external computer-based data collection systems.

Electrical

Monitor: 12 VDC Power, max. 3 A, by external power supply; Type BF, Applied Part, Defibrillation Proof.

Power supply for hospital use: Class II FE (with functional earth), Electrical Safety (IEC 60601-1)

Power supply for home use: Class II (without functional earth), Electrical Safety (IEC 60601-1)

Internal battery type: rechargeable, sealed Li Ion Battery

Battery Capacity (new, fully-charged): Normal Operation: up to 4 hours; Sleep Mode: up to 7 hours

Charging Time: approximately 4 hours

Environmental

Transport / storage temperature: 0 – 50°C (32 – 122°F)

Transport / storage humidity: 10 – 90% non-condensing

Operating temperature: 5 – 40°C (41 – 104°F)

Operating humidity: 15 – 90% non-condensing

Operating altitude: -400 – 5000 m (-1300 – 16404 ft)

Accuracy: ± 3 mmHg (0.4 kPa)

Built-in barometer: Range: 350 – 820 mmHg (47 – 109 kPa); Accuracy: ± 3 mmHg (0.4 kPa)

Compliance

IEC 60601-1, ANSI/AAMI ES60601-1, CAN/CSA-C22.2 No. 60601-1, IEC 60601-1-6, IEC 60601-1-8, IEC 60601-2-23, ISO 80601-2-61, 60601-1-11, 60601-1-2

Regulatory Approval

The tCOM+ monitor is FDA-cleared and fulfills the General Safety and Performance Requirements stipulated by the EU Medical Device Regulation (EU) 2017/745.