



ASPIRATOR

Technical features

- The electric aspirator can suction the pus, blood, sputum and other glue liquid effectively, and meet the most requirements for operating room, ward department, outdoor rescue, etc.
- This product is portable, and use the durable plastic shell with 12V5AH battery inside, and can be used with vehicle powersupply DC12V or charge by the city electric AC220V.



Technical parameter

- Extreme pressure bearing value: 650mmHg
- Pressure bearing adjustable scope: 150~650mmHg
- Air extracting rate: 25L/min
- Noise: $\leq 60\text{dB(A)}$
- Liquid storage bottle capacity: 500ml
- Power supply: AC220V $\pm 10\%$, 50Hz or DC12V
- Self-weight: 5KG
- Size: 35 \times 28 \times 16cm

Accessories

