

Product Specifications

Suction Systems with Regulators

Double suction - 2 x 2L polycarbonate canisters which graduate in 100ml increments, with a regulator, cage and Heyer connection

Single suction - 2L canister with regulator, cage and brass Heyer connection

Single suction - 1L canister with cage and tubing and Heyer connection

Regulated vacuum controllable vacuum power

Robust design with high pressure regulator

Low pressure regulator available on request

Overflow protected

Inline filter optional

Can be supplied for direct mount, gabler/accessory mount or on portable trolley

Used in conjunction with disposable liners

On/off switch

Operating pressure 0-100kpa/mmHg

Colour coded vacuum hose - 2m length

Flow rate >25l/m

Hydrophobic filter included

Suction inlet and connectors and chrome plated brass

Gauge is marked in kPa and mmHg

Gauge is fully removable for easy maintenance

Easy to service

Pipeline Suction Units Single / Double



Suction Liners With Filter / Without Filter



1ℓ

2ℓ

2ℓ

Suction Canisters 1ℓ / 2ℓ

