

# Global Markets



**Amoul** Ambulanc (Shenzhen) Tech. Co., Ltd.



Building #8, Skyworth Innovation Valley Phase II, TangTou Road No.1, ShiYan, Bao'An District, ShenZhen, China  
+86 755-26072286  
www.amoulmed.com  
info@amoulmed.com

Amoulmed ambulshenzhen.Tech Amoulmed Amoulmed Amoul

© Copyright 2024 Amoulanc(Shenzhen) Technology Co., Ltd  
ISBN 9787552607228 Printed in China

**Amoul**<sup>®</sup>

# External Defibrillator

(Fully automated/semi-automated)



www.amoulmed.com

# AMOUL® AED — i-series

New Pioneer in Public Emergency Rescue

## i7/i9

	Screen	Types
i7	Graphic Indication	Semi-Automated
i7 Plus	Graphic Indication	Automated
i9	7" Color Screen	Semi-Automated
i9 Plus	7" Color Screen	Automated



### Higher in Energy

Adopting biphasic exponential truncated wave  
The defibrillation shocks are delivered in a progressive energy sequence to 360J



### Quick in Response

Simultaneous performing of heart rate analysis and charging, ensuring efficient rescue



### Easy in Use

One-step Power ON



Clear display with 7" HD screen



One-key switch patient type



Fast in language selection

### Elegant Design



Strong Appearance



High Protection Class IP65



Large Storage Capacity

# Safe Management and Comprehensive Protection

AED Intelligent Network Management Platform



Remote device management



Remote maintenance of consumables



Remote data viewing

## AED Cabinet




Wall Cabinet (Metal)



Wall Cabinet (Plastic)



Vertical Cabinet



## First-aid Training Solutions



AED Trainer



Blood Flow Visualization Manikin



Perfect CPR Training Solution