

yuwell

BREATHCARE III AUTO CPAP

YH-680A/B

Color ● Black ● White

Mode CPAP/Auto CPAP

Noise 27±2 (dBA)

Pressure 4-20 (cmH₂O)



Professional Care

- Intelligent Algorithm
- Multiple-events Monitoring
- Ultra-quiet Noise
- Efficient Humidifier
- Central Apnea Detection
- Mask-fit Test

Auto-all

- Auto Ramp
- Auto Heated Tube
- Auto Light
- Smart Start/Stop
- Altitude Compensation
- Leak Compensation

Wireless Data System

- Cellular/WiFi Module
- Free Data-cloud Access
- OTA
- ☁ BreathCare Web
- 📺 BreathCare Station
- 📱 BreathCare+

Packaging parameters

Model
YH-680A/B

QTN
(PCS/CTN)
1

G.W./N.W.
(kgs)
3.7/2.4

Size
(cm)
41x31x18

Ecomed
Medical (Pty) Ltd

Ecomed Medical (Pty) Ltd
www.ecomed.co.za

25 Clew Street, Monument Ext. 1, Krugersdorp, 1739
+27 11 955 5710 | cpap@ecommed.co.za