

BREATHWEAR MASK SERIES

YF-02



Fit Comfort Convenience Light & Soft

- Type: Full Face Mask
- Therapy Pressure Range :4-30 cmH2O
- Tube Connection : \varnothing 22mm
- Flow Resistance(At 50 L/min): \leq 0.5cmH2O
- Flow Resistance(At 100 L/min): \leq 1.5cmH2O
- Expiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min): \leq 2.0cmH2O/(L/S)
- Inspiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min): \leq 2.0cmH2O/(L/S)
- Anti Asphyxia Valve Open Pressure: \leq 2.0cmH2O/(L/S)
- Anti Asphyxia Valve Closed Pressure: \leq 2.0cmH2O/(L/S)



PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
YF-02	10	2/1.6	35×28×28