

# POSITIVE AIRWAY PRESSURE UNITS

## YH-350

- Mode: CPAP
- Pressure: 4-20cmH<sub>2</sub>O
- Noise: 28dBA±2dBA
- Mask (Optional)
- SD card (Optional)
- Smart Start/Stop
- Altitude Adjustment
- Leakage Compensation
- Compliance Report
- Data Management
- Reminder Function
- Ramp Function
- Efficient Humidifier



## PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
YH-350	1	2.9/1.5	41.5×30×16