

# BREATHWEAR MASK SERIES

## YP-01



Fit



Comfort



Convenience



Light & Soft

- Type: Pillow mask
- Therapy Pressure Range :4-20 cmH2O
- Tube Connection :  $\varnothing$ 22mm
- Flow Resistance(At 50 L/min):  $\leq$ 0.5cmH2O
- Flow Resistance(At 100 L/min):  $\leq$ 1.5cmH2O
- Expiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min):  $\leq$ 2.0cmH2O/(L/S)
- Inspiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min):  $\leq$ 2.0cmH2O/(L/S)
- Anti Asphyxia Valve Open Pressure:  $\leq$  2.0cmH2O/(L/S)
- Anti Asphyxia Valve Closed Pressure:  $\leq$  2.0cmH2O/(L/S)



## PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
YP-01	10	3/2.82	43×21.5×16.8