

Regulator - Bullnose



**Dial Click Stop Constant Flow with pressure gauge to register cylinder pressure
Oxygen delivery through Click-Stop
Hand Dial Marked 0-25l/min
Oxygen Cylinder Pressure Gauge
Pressure Dial protected by a rubber ring
Outer Body manufactured from Aluminium
Inner Body Manufactured from Brass
Minimum Input Operating Pressure 200 bars
Output Pressure at Heier: 3 - 6 Bar
Complete Male/Female Connection Included**

2 Year Warranty & 3x O-Rings included