

# Ecomed Medical (Pty) Ltd

## Oxygen Flowmeters

Precision engineered pressure compensated design accurately measures flow from 0 to 15  $\ell$ pm  
Easy to read flow tube and unbreakable transparent polycarbonate outer cover for strength and visibility  
Chromed brass body or anodized aluminum body for durability  
Stainless steel needle valve for precise control.  
Calibrated at 50 psi inlet pressure  
Outlet flow control valve 180°  
Inlet and outlet thread size are available for International standard

### Ultra Flowmeter Single Ball & Tube

- Ultra Single Ball & Tube
- Inlet Pressure:  $\leq 400$  Bar
- Increments:  
1,2,3,4,5,6,7,8,9,10,11,12,13,14,15  $\ell$ pm



### Ultra Flowmeter Double Ball & Tube

- Ultra Single Ball & Tube
- Inlet Pressure:  $\leq 400$  Bar
- Increments:  
1,2,3,4,5,6,7,8,9,10,11,12,13,14,15  $\ell$ pm

