

BREATHWEAR MASK SERIES

YN-03



Fit Comfort Convenience Light & Soft

- Type: Nasal Mask
- Therapy Pressure Range :4-30 cmH2O
- Tube Connection : φ 22mm
- Flow Resistance(At 50 L/min): \leq 0.5cmH2O
- Flow Resistance(At 100 L/min): \leq 1.5cmH2O



PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
YN-03	10	1.7/1.3	35×28×28