



POSITIVE AIRWAY PRESSURE UNITS

YH825

- Mobile App Connectivity
- Mode: CPAP, S, T, ST, VGPS
- Pressure: 4 - 25cmH₂O
- Effective Humidifier
- High Pressure Reminder
- Leak Reminder
- Asphyxiation Reminder
- Trigger Adjustment
- Cycle Adjustment
- Slope Adjustment
- Tidal Volume Adjustment
- Altitude Adjustment
- Leak Compensation
- Data Management
- Central Apnea Detection
- Smart Humidification

