



Assisted Reproductive Medicine

Voluson SWIFT This Changes Everything


Your success is defined by life. Prospective parents seek out your expertise in order to conceive. That's why the Voluson™ SWIFT ultrasound system delivers exceptional imaging and simple yet innovative tools to help you and your patients achieve their dream of a successful pregnancy.




Why Can't We Conceive? Infertility Assessment

Effortless Imaging

Exceptional imaging is essential when investigating causes of infertility and Voluson Core Architecture enables exceptional imaging in any mode.




Normal Uterus




Power Doppler Endometrioma

3D Simplified

Adding 3D to your practice can provide additional information in the assessment of infertility and now it is simple with Voluson SWIFT.




2D to 3D coronal uterus in 3 steps with Uterine Trace



Ovarian mass displayed using Tomographic Ultrasound Imaging (TUI)




HyCoSy[®] (Hysterosalpingo Contrast Sonography) demonstrating patent fallopian tubes




SonoAVC[™] antral^{2.0} Sonography-based Automated Volume Count automates ovarian reserve assessment

Improving Efficiency in Follicle Monitoring

Follicle monitoring is one of the most critical aspects of fertility management and can be the most time consuming. Automate follicle counts in 2D or 3D with Auto Caliper and SonoAVC^{follicle} helping increase efficiency and consistency.



Auto Caliper – Tap follicle to measure



SonoAVC^{follicle} automatically calculates follicle volume and dimensions

Safe, Confident Oocyte Retrieval and Embryo Transfer

Voluson SWIFT delivers in every way during procedures

- Clear visualization during oocyte retrieval and embryo transfer
- Optimized user and patient comfort with the IC9B-RS bent endocavity probe
- Full Screen Mode for even easier visualization
- Simplified design to aid in cleaning and disinfection



Oocyte Retrieval

Confirming a Dream

Identify an early pregnancy with ease. With exceptional imaging easily demonstrate fetal heart activity to parents. As the pregnancy progresses adding 3D can add a more realistic view.



8 Week 5 Day Fetus with Yolk Sac



HDlive[™] 11-week Fetus

*HyCoSy is not available in all countries.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

©2023 GE HealthCare.

Voluson, SonoAVC and HDlive are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license. March 2023 JB22802XX