

Ecomed

Medical (Pty) Ltd

Oxygen Cylinders

Product Specifications

Diameter: 11cm

Height: 61cm

Weight (Empty): 3.54kg

Oxygen Litres: 420ℓ

Water Litres: 3.2ℓ

Flow @ 2L/min: 3 Hours 30 Minutes

Pin Index Oxygen Cylinder Size D / Portable



Specialists in Home Oxygen, CPAP, Hospital, Trauma and Pain Management Equipment

Tel: 011 955 5710 Fax: 011 955 6316 Email: sales@ecomед.co.za