

TRIPLE YOUR OXYGEN CYLINDER LIFE

EASYPULSE5 Oxygen Conserving Regulator

- Extends Cylinder Life by Three to Six Times
*Depending on Patients' Set Flow Rates
- Reduces Oxygen Refill Costs
- Light Weight & Portable
- Easy Cylinder Attachment
- Ideal for **COPD & Long-Term Oxygen Therapy**



EASYPULSE5 Oxygen Conserving Regulator

Model: 198705

Flow Setting	Pulse Volume @ Breath Rate of:				Minute Volume
	15	20	25	30	
1	17ml	13ml	10ml	9ml	260ml/min
2	28ml	21ml	17ml	14ml	420ml/min
3	37ml	28ml	22ml	19ml	560ml/min
4	47ml	35ml	28ml	23ml	700ml/min
5	57ml	43ml	34ml	29ml	860ml/min
2lpm	N/A				2000ml/min (2lpm)

Minute Volume Consistent @12 = 38 BPM

Cylinder Duration		Flow Setting				
Cylinder Size	O2 Capacity	1	2	3	4	5
Size D - Portable	420 liters	25hrs	16hrs	12hrs	10hrs	8hrs
Size F - Back-Up	1180 liters	70hrs	45hrs	34hrs	28hrs	22hrs

The **EasyPulse5** oxygen conserving regulator delivers oxygen on inhalation, helping extend cylinder duration and reducing refill costs.