



Homecare Product Catalogue

Ecomed Medical's Homecare Catalogue brings essential health products right to your home.

Shop at one of our 23 branches nationwide or through our online store where we make it easy to access the right products and support whenever you need them.

www.ecomed.co.za

Shop Online





This division of Ecomed was established in the year 2000 with a focus on supplying comprehensive oxygen and CPAP solutions to all medical schemes across South Africa. Over the years, we have grown to serve more than 10,000 patients who currently benefit from our high-quality oxygen concentrators, ensuring reliable and effective therapy in the comfort of their homes.

We have a dedicated team of homecare advisors who visit patients personally, providing guidance, support, and training on the correct use of their equipment. All equipment is serviced and maintained through our two fully equipped technical centers, ensuring reliability and safety at all times.

In addition to oxygen and CPAP equipment, we offer flexible rental options for oxygen cylinders, pulse oximeters, mobile suction units, and humidifiers, allowing patients to access the right equipment according to their individual needs. With our extensive experience and nationwide reach, Ecomed remains a trusted provider of oxygen and CPAP solutions to all medical aids in South Africa, committed to improving respiratory health and enhancing patient care.

Products

Offering a wide range of quality homecare products.

OXYGEN THERAPY

Stationary Oxygen Concentrators Page 4-5

Portable Oxygen Concentrators Page 6-7

Oxygen Cylinders Page 8

HOMECARE EQUIPMENT

OPEP Device Page 10

Blood Pressure Monitors Page 11

Nebulisers Page 12

Suction Units Page 13

Pulse Oximeters Page 14



OXYGEN THERAPY

Stationary Oxygen Concentrators

A patient would use a stationary oxygen concentrator to receive a continuous, reliable supply of oxygen at home, supporting long-term therapy and stable oxygen levels while resting, sleeping, or carrying out daily activities indoors.

Portable Oxygen Concentrators

A patient would use a portable oxygen concentrator to maintain safe oxygen levels while staying mobile and independent, allowing them to move freely, travel, and carry out daily activities without relying on heavy oxygen cylinders.

Oxygen Cylinders

A patient would use oxygen cylinders to receive a dependable supply of medical oxygen when a concentrator is not available, providing immediate oxygen therapy at home, in transit, or during emergencies.



5L Oxygen Concentrator

Ecomed Medical's oxygen concentrators deliver safe, uninterrupted oxygen therapy in the comfort of your home.

Easy to use, energy-efficient, and supported by our nationwide service network, they ensure you or your loved ones breathe easier every day.



Yuwell 5L Oxygen Concentrator



High Purity



Failure Detection



Ultra-Silent



Castors

Provides high-purity oxygen at continuous flow up to 5L/min, featuring safety alarms for low oxygen, power failure, and pressure issues, quiet operation under 52dB for home comfort, easy-to-use LCD display with simple controls and built-in humidifier, and durable, reliable construction for long-term dependable use.

Buy Now



2 Year Warranty


Accessories



Humidifier Bottle



Nasal Cannula



Extension Tubing



Oxygen Mask

10L Oxygen Concentrator

Ecomed Medical's oxygen concentrators deliver safe, uninterrupted oxygen therapy in the comfort of your home.

Easy to use, energy-efficient, and supported by our nationwide service network, they ensure you or your loved ones breathe easier every day.



Yuwell 10L Oxygen Concentrator



High Purity



Failure Detection



Ultra-Silent



Castors

Delivers powerful performance with up to 10L/min continuous oxygen flow and consistent purity of approximately 87%–95.5%, featuring smart safety alarms for power issues, low oxygen, and pressure problems, quiet and easy-to-use controls for comfortable daily home use, and a durable, convenient design with filter reminders and caster wheels for long-term reliability and mobility.

Buy Now



2 Year Warranty


Accessories



Humidifier Bottle



High Flow Nasal Cannula



Extension Tubing



Oxygen Mask

Portable Oxygen Concentrators

Ecomed Medical's portable oxygen concentrators provide reliable, continuous oxygen therapy wherever you are.

Compact, user-friendly, and energy-efficient, they make daily breathing easier—backed by our nationwide service network for support and peace of mind.



Kingon P2 Portable Oxygen Concentrator



High Purity



Safety Alarms



Light Weight



Pulse Dose



Triple Filtration

This ultra-light and portable device, weighing under 2kg and FAA approved for safe air travel, delivers efficient, breath-synchronized pulse-dose oxygen, offers long battery life of up to 4 hours with both AC and DC power options, and features a smart, easy-to-use LCD display with built-in safety alarms for convenient and reliable use on the go.


Buy Now



3 Year Warranty

Exclude: Battery, Sieve Beds, Power Supply

Accessories




Backpack



Extra Battery



Battery Charger



AC Charger



DC Charger

Portable Oxygen Concentrators

Ecomed Medical's portable oxygen concentrators provide reliable, continuous oxygen therapy wherever you are.

Compact, user-friendly, and energy-efficient, they make daily breathing easier—backed by our nationwide service network for support and peace of mind.



Inogen G5 Portable Oxygen Concentrator



High Purity



Safety Alarms



Light Weight



Pulse Dose



Ultra-Silent

This ultra-light and portable device, weighing approximately 2.2kg for easy travel, offers long battery life of up to 13hours with an extended battery, features adjustable pulse-dose settings from 1 to 6 for personalized oxygen therapy, is FAA-approved for safe air travel, and combines quiet operation with user-friendly, easy-to-use controls for convenient, reliable use anytime and anywhere.

Buy Now



3 Year Warranty

Exclude: Battery, Sieve Beds, Power Supply

Accessories



Carry Bag

Oxygen Cylinders

Ecomed Medical's oxygen cylinders provide safe, reliable, and continuous oxygen therapy wherever it's needed.

Available in a range of sizes, they are easy to handle, transport, and use, ensuring uninterrupted breathing support.

Cylinders can be conveniently refilled at any of our 23 branches nationwide, and with our full range of compliant accessories.



Portable 424L Oxygen Cylinder Kit



This device measures 11cm in diameter and 61cm in height, weighs 3.54kg when empty, holds 420 liters of oxygen and provides a continuous oxygen flow of 2L/min for approximately 3 hours and 30 minutes.

Buy Now



Back-Up 1400L Oxygen Cylinder Kit



This device has a diameter of 14cm and a height of 82cm, weighs 7.86kg when empty, holds 1,120 liters of oxygen and provides a continuous oxygen flow of 2L/min for up to 9 hours.

Buy Now



HEMOCARE EQUIPMENT

OPEP Device

A patient would use an OPEP device to help loosen and clear mucus from the airways, improving breathing, reducing congestion, and supporting better lung function.

Blood Pressure Monitors

A patient would use a blood pressure monitor to regularly track their blood pressure at home, helping detect changes early, manage chronic conditions, and support effective treatment.

Nebuliser

A patient would use a nebuliser to deliver medication directly to the lungs in a mist form, making it easier to treat respiratory conditions like asthma, COPD, or bronchitis and improve breathing.

Suction Unit

A patient would use a home suction unit to safely remove mucus and secretions from the airway, helping maintain clear breathing and prevent respiratory complications at home.

Pulse Oximeter

A patient would use a pulse oximeter to quickly and non-invasively monitor their blood oxygen levels and heart rate, helping detect low oxygen early and manage respiratory or cardiovascular conditions effectively.



OPEP Device

The Trudell Aerobika is a handheld oscillating PEP device designed to help clear mucus from the airways.

It uses gentle vibrations and positive pressure during exhalation to loosen secretions, making breathing easier and supporting effective airway clearance for patients with respiratory conditions.



Trudell Aerobika OPEP Device



Drug Free



Adjustable Settings



No Side Effects



Reusable

This clinically proven airway clearance device uses oscillating PEP therapy to help loosen and mobilize mucus for easier expectoration, features adjustable resistance levels to allow personalized therapy tailored to different patient needs and lung capacities, is drug-free and easy to use with a simple handheld design that supports effective daily therapy at home or in clinical settings, works in any position—sitting, standing, or lying down—to improve patient compliance, and is built to be durable and easy to clean, with simple disassembly for proper hygiene and repeated use.

Buy Now





How the *Aerobika*® device works.


The *Aerobika*® device has been clinically proven to improve ventilation, lung function and quality of life as well as reduce COPD flare-ups.¹




BEFORE USE: Airways inflamed, collapsed and plugged with mucus



AFTER USE: As you exhale through the device, the airways are opened



As the airways open, vibrations help thin and loosen mucus which can



This may help your inhaled medication work better²

Blood Pressure Monitors

Ecomed Medical's blood pressure monitors provide accurate, reliable readings anytime, anywhere. Compact, easy-to-use, and clinically validated, they help you keep track of your health—backed by our nationwide service network for support and peace of mind.



Yuwell Upper Arm Blood Pressure Monitor



One Touch Operation



Memory Storage



Cuff Detection



Irregular Detection



Audio Feedback

The Yuwell Upper Arm Blood Pressure Monitor provides accurate and reliable measurements of blood pressure and pulse, featuring an easy-to-read LCD display, one-touch automatic operation, memory storage for multiple readings, irregular heartbeat detection, and a comfortable adjustable cuff, making it ideal for convenient daily home monitoring and effective health management.

Buy Now



1 Year Warranty

Yuwell Wrist Blood Pressure Monitor



One Touch Operation



Memory Storage



Cuff Detection



Irregular Detection

The Yuwell Wrist Blood Pressure Monitor provides fast and accurate blood pressure and pulse measurements, featuring a large, easy-to-read LCD display, one-touch automatic operation, memory storage for multiple readings, irregular heartbeat detection, adjustable wrist cuff for comfort, and a compact, portable design, making it ideal for convenient daily monitoring at home or on the go.

Buy Now



1 Year Warranty

Nebulisers

Ecomed Medical's nebulisers provide safe, effective, and reliable respiratory therapy in the comfort of your home or wherever you need it.

Designed for ease of use, they are compact, lightweight, and quiet, making treatment simple for patients of all ages.



Yuwell Mesh Nebuliser



Dual Power



Ultra-Silent



Minimal Waste



Multiple Accessories

Featuring advanced mesh technology for fast and efficient medication delivery, this nebuliser is portable and lightweight for use anywhere, operates quietly for comfort even during sleep, minimises medication waste to maximise drug use, and offers flexible power options with battery or USB operation.

Buy Now



1 Year Warranty

Yuwell Piston Nebuliser



Deep Lung Delivery



Ultra-Silent



Minimal Waste



Oil Free

This nebuliser delivers a fine mist that reaches deep into the lungs for effective treatment, features a reliable, oil-free compressor for consistent airflow, is compact, lightweight, and easy to use, operates quietly for comfort across all ages, and minimises medication waste to maximise drug use.

Buy Now



1 Year Warranty

Suction Units

Ecomed Medical's suction units deliver safe, efficient, and reliable airway clearance therapy wherever it's needed —whether at home, in a healthcare setting, or on the go.

Designed with powerful yet gentle suction, they help remove mucus and secretions effectively, promoting easier breathing and reducing the risk of respiratory complications.



Yuwell Battery Operated Suction



- Strong Suction
- Dual Power
- Safety Features
- 1000ml Jar
- Ultra-Silent

Featuring dual power operation on battery or AC/DC, this suction unit delivers powerful suction to efficiently remove mucus and secretions, is portable and durable with a lightweight design and built-in handle, incorporates hygienic overflow protection and bacterial filters, and operates quietly with simple, easy-to-use controls for comfortable use.

Buy Now



2 Year Warranty

Yuwell Electric Suction



- Strong Suction
- Oil Free
- Safety Features
- 1000ml Jar
- Ultra-Silent

This suction unit delivers powerful, adjustable suction for effective mucus removal, features a maintenance-free, oil-free pump for reliable and easy care, includes safety features such as overflow protection and bacterial filters, offers user-friendly front panel controls, and combines quiet operation with a durable, long-lasting design.

Buy Now



2 Year Warranty

Pulse Oximeter

Ecomed Medical's pulse oximeters provide fast, accurate, and reliable monitoring of blood oxygen levels and heart rate — anywhere, at home or on the go.

Compact, lightweight, and easy to use, they allow patients and caregivers to track vital signs with confidence, helping detect changes early and supporting better overall health management.



Biolight M70 Finger Tip Pulse Oximeter



One Step Operation



Auto Power Off



Battery Powered



Accurate

This pulse oximeter provides fast and accurate SpO₂ and pulse readings, features a clear OLED display that's easy to read even in low light, includes auto power-off to save battery life, is portable and lightweight for home or travel use, and is built with a durable, shock-proof design for long-lasting reliability.

Buy Now



1 Year Warranty

Biolight M800 Handheld Pulse Oximeter



Baby



Dual Monitoring



Visual & Audio Alerts



Memory Storage

This versatile monitoring device measures SpO₂, pulse, ECG, and respiratory rate, features a large, easy-to-read colour OLED display, supports both spot checks and continuous monitoring, offers long battery life of up to 16 hours, and is portable and lightweight for convenient use at home, in clinics, or on the go.

View Brochure



1 Year Warranty

GAUTENG

Head Office - Krugersdorp
 25 Clew Street, Monument Ext.1, 1745
 Tel: 011 955 5710 / On-Call: 082 729 0749
 Office Hours: Mon- Fri 08:00 - 16:30

Master Warehouse
 Unit 5, Block B, Factoria Industrial Estate,
 25 Louis Friedman Street, Factoria, 1739
 Mon - Fri 08:00 - 16:30, Sat - Sun 08:00 - 12:00

Homecare Technical
 Unit 9, Block C, Factoria Industrial Estate,
 25 Louis Friedman Street, Factoria, 1739
 Mon- Fri 08:00 - 16:30

East Rand
 Unit 23, Oakhurst Office Park, 5 Grippen Rd, 1459
 Tel: 011 027 4259
 Mon - Fri 08:00 - 16:30, Sat 08:00 - 13:00

Pretoria Moot
 271 15th Avenue, Rietfontein, 0084
 Tel: 012 993 0493
 Mon - Fri 08:00 - 16:30, Sat 08:00 - 13:00

Pretoria East
 Unit 35, Scientia Techno Park,
 2 Meiring Naudé Rd, Lynnwood, 0157
 Tel: 012 885 0890
 Mon- Fri 08:00 - 16:30

Vaal
 137 General Hertzog Rd, Three Rivers,
 Vereeniging, 1935
 Tel: 016 013 0090
 Mon- Fri 08:00 - 16:30

Centurion
 13 Esdoring Nook, Highveld Technopark,
 Centurion, 0157
 Tel: 086 111 4028
 Mon- Fri 08:00 - 16:30

FREE STATE

Bloemfontein
 9 Nig Mientjie Street, Pellissier, 9301
 Tel: 051 422 4118
 Mon- Fri 08:00 - 16:30

Welkom
 144 Constantia Street, Dagbreek, 9459
 Tel: 057 030 0130
 Mon- Fri 08:00 - 16:30

Bethlehem
 43 De Leeuw St., Hospitaalheuwel, 9701
 Tel: 082 593 8068
 Mon- Fri 08:00 - 13:00

MPUMALANGA

Nelspruit
 Unit 10 Central Park, Suikerriet Street, 1200
 Tel: 010 446 9832
 Mon- Fri 08:00 - 16:30

Witbank
 1 Maria Street, Fransville, 1035
 Tel: 013 590 0559
 Mon- Fri 08:00 - 16:30

KWA-ZULU NATAL

Durban
 450A Che Guevara, Glenwood, 4083
 Tel: 031 202 4214
 Mon- Fri 08:00 - 16:30

Dundee
 84 Victoria Street, 3000
 Tel: 082 491 3515
 Mon- Fri 08:00 - 16:30

Pietermaritzburg
 296 Bulwer St, CBD, 3201
 Tel: 067 769 9057
 Mon- Fri 08:00 - 13:00

WESTERN CAPE

Milnerton
 Unit 36 Frazzitta Business Park, Freedom Way,
 Milnerton, Cape Town, 7441
 Tel: 021 020 1750
 Mon- Fri 08:00 - 16:30

Muizenburg
 4 Tokai Close, Marina de Gama, 7945
 Tel: 021 788 3163
 Mon- Fri 08:00 - 16:30

George
 22 Waboom St, Glen Barrie, George, 6529
 Tel: 044 124 0160
 Mon- Fri 08:00 - 16:30

NORTH WEST

Klerksdorp
 Unit 7, Stargate Office Park,
 20 Buffeldoorn Avenue, Wilkopies, 2571
 Tel: 010 446 9833
 Mon- Fri 08:00 - 16:30

Rustenburg
 Unit 4, 250 Beyers Naude Dr,
 Rustenburg, 2999
 Tel: 014 112 0043
 Mon- Fri 08:00 - 16:30

EASTERN CAPE

Port Elizabeth
 29 Starling Crescent, Greenshields Park,
 Gqeberha, 6070
 Tel: 086 111 1718
 Mon- Fri 08:00 - 16:30


East London
 75 Devereux Avenue, Vincent, 5241
 Tel: 010 446 9834
 Mon- Fri 08:00 - 16:30

LIMPOPO

Polokwane
 56 Bendor Drive, Bendor, 0699
 Tel: 010 446 9831
 Mon- Fri 08:00 - 16:30

NORTHERN CAPE

Kimberley
 14A Schmidtdrift Rd, Verwoerdpark, 8301
 Tel: 072 141 2335
 Mon- Fri 08:00 - 16:30



Ecomed

Medical (Pty) Ltd



Through our 23 branches nation wide, we ensure patients have continuous, uninterrupted therapy.



+27 11 955 5710

We offer a 24-hour helpline, ensuring a patient's needs are supported at any time, day or night.

Visit our website



www.ecomed.co.za