



Homecare Product Catalogue

Ecomed Medical's Homecare Catalogue brings essential health products right to your home.

Shop at on of our 24 branches nationwide or through our online store where we make it easy to access the right products and support whenever you need them.

www.ecomed.co.za

Shop Online





This division of Ecomed was established in the year 2000 with a focus on supplying comprehensive oxygen and CPAP solutions to all medical schemes across South Africa. Over the years, we have grown to serve more than 10,000 patients who currently benefit from our high-quality oxygen concentrators, ensuring reliable and effective therapy in the comfort of their homes.

We have a dedicated team of homecare advisors who visit patients personally, providing guidance, support, and training on the correct use of their equipment. All equipment is serviced and maintained through our two fully equipped technical centers, ensuring reliability and safety at all times.

In addition to oxygen and CPAP equipment, we offer flexible rental options for oxygen cylinders, pulse oximeters, mobile suction units, and humidifiers, allowing patients to access the right equipment according to their individual needs. With our extensive experience and nationwide reach, Ecomed remains a trusted provider of oxygen and CPAP solutions to all medical aids in South Africa, committed to improving respiratory health and enhancing patient care.

Products

Offering a wide range of quality homecare products.

OXYGEN THERAPY

Stationary Oxygen Concentrators Page 4-5

Portable Oxygen Concentrators Page 6-7

Oxygen Cylinders Page 8

Regulators Page 9

HEMECARE EQUIPMENT

OPEP Device Page 11

Spacer Device Page 12-13

Blood Pressure Monitor Page 14

Nebulisers Page 15-16

Suction Units Page 17

Pulse Oximeters Page 18

OXYGEN THERAPY

Stationary Oxygen Concentrators

Stationary oxygen concentrators are designed to deliver a continuous, reliable supply of oxygen at home, supporting long-term therapy and stable oxygen levels while resting, sleeping, or carrying out daily activities indoors.

Portable Oxygen Concentrators

A portable oxygen concentrator maintains safe oxygen levels while staying mobile and independent, allowing them to move freely, travel, and carry out daily activities without relying on heavy oxygen cylinders.

Oxygen Cylinders

Oxygen cylinders deliver a dependable supply of medical oxygen when a concentrator is not available, providing immediate oxygen therapy at home, in transit, or during emergencies.

Oxygen Regulators

An oxygen regulator is a device that attaches directly to an oxygen cylinder and is used to control and reduce the pressure of the oxygen being delivered by the cylinder to a safe, usable level.



5L Oxygen Concentrator

Ecomed Medical's oxygen concentrators deliver safe, uninterrupted oxygen therapy in the comfort of your home.

Easy to use, energy-efficient, and supported by our nationwide service network, they ensure you or your loved ones breathe easier every day.



Yuwell 5L Oxygen Concentrator



High Purity



Failure Detection



Ultra-Silent



Humidifier

Provides high-purity oxygen at continuous flow up to 5L/min, featuring safety alarms for low oxygen, power failure, and pressure issues, quiet operation under 52dB for home comfort, easy-to-use LCD display with simple controls and built-in humidifier, and durable, reliable construction for long-term dependable use.

Buy Now



2 Year Warranty


Accessories



Humidifier Bottle



Nasal Cannula



Extension Tubing



Oxygen Mask

10L Oxygen Concentrator

Ecomed Medical's oxygen concentrators deliver safe, uninterrupted oxygen therapy in the comfort of your home.

Easy to use, energy-efficient, and supported by our nationwide service network, they ensure you or your loved ones breathe easier every day.



Yuwell 10L Oxygen Concentrator



High Purity



Failure Detection



Ultra-Silent



Humidifier

Delivers powerful performance with up to 10L/min continuous oxygen flow and consistent purity of approximately 87%—95.5%, featuring smart safety alarms for power issues, low oxygen, and pressure problems, quiet and easy-to-use controls for comfortable daily home use, and a durable, convenient design with filter reminders and caster wheels for long-term reliability and mobility.

Buy Now



2 Year Warranty


Accessories



Humidifier Bottle



High Flow Nasal Cannula



Extension Tubing



Oxygen Mask

Portable Oxygen Concentrators

Ecomed Medical's portable oxygen concentrators provide reliable, continuous oxygen therapy wherever you are.

Compact, user-friendly, and energy-efficient, they make daily breathing easier—backed by our nationwide service network for support and peace of mind.



Kingon P2 Portable Oxygen Concentrator



High Purity



Safety Alarms



Light Weight



Pulse Dose



Triple Filtration

This ultra-light and portable device, weighing under 2kg and FAA approved for safe air travel, delivers efficient, breath-synchronized pulse-dose oxygen, offers long battery life of up to 4 hours with both AC and DC power options, and features a smart, easy-to-use LCD display with built-in safety alarms for convenient and reliable use on the go.


Buy Now



3 Year Warranty

Exclude: Battery, Sieve Beds, Power Supply

Accessories




Backpack



Extra Battery



Battery Charger



AC Charger



DC Charger

Portable Oxygen Concentrators

Ecomed Medical's portable oxygen concentrators provide reliable, continuous oxygen therapy wherever you are.

Compact, user-friendly, and energy-efficient, they make daily breathing easier—backed by our nationwide service network for support and peace of mind.



Inogen Rove 6 Portable Oxygen Concentrator



High Purity



Safety Alarms



Light Weight



Pulse Dose



Ultra-Silent

The Inogen Rove 6 Portable Oxygen Concentrator is a lightweight (± 2.2 kg), travel-friendly device with powerful oxygen output, offering up to ± 6.5 hours of battery life with the standard battery and up to ± 12.75 hours with the extended battery, along with quiet operation and smart breath-triggered delivery—making it a reliable, low-maintenance solution for oxygen therapy on the go.

Buy Now



3 Year Warranty

Exclude: Battery, Sieve Beds, Power Supply

Accessories



Carry Bag



Extended Battery



Battery Charger



AC Charger



DC Charger

Oxygen Cylinders

Ecomed Medical's oxygen cylinders provide safe, reliable, and continuous oxygen therapy wherever it's needed.

Available in a range of sizes, they are easy to handle, transport, and use, ensuring uninterrupted breathing support.

Cylinders can be conveniently refilled at any of our 24 branches nationwide, and with our full range of compliant accessories.



Portable (D) Oxygen Cylinder Kit



The cylinder measures 11cm in diameter and 61cm in high, weighs 3.54kg when empty, holds 420 liters of oxygen and provides a continuous oxygen flow of 2L/min for approximately 3 hours and 30 minutes.

Buy Now



Back-Up (F) Oxygen Cylinder Kit



The cylinder has a diameter of 14cm and a height of 82cm, weighs 7.86kg when empty, holds 1,120 liters of oxygen and provides a continuous oxygen flow of 2L/min for up to 9 hours.

Buy Now



Oxygen Regulators

Ecomed Medical's oxygen cylinders are fitted with high-quality pin index regulators and conserving regulators, ensuring safe, accurate, and efficient oxygen delivery at all times.

Designed for reliability and ease of use, these systems help optimise oxygen usage while maintaining consistent therapy, whether at home or on the go.



Pin Index Regulator



The Ultra Pin Index regulator operates with an inlet pressure of 200 Bar and delivers an outlet pressure of 3–6 Bar, offering adjustable flow rate increments of 1, 2, 4, 6, 8, 10, 12, 15, 18, 20, and 25 LPM.

Buy Now



Conserving Regulator



The EasyPulse conserving regulator is designed to extend cylinder life by delivering oxygen only when the patient inhales, helping you get more usage time from each cylinder while maintaining effective oxygen therapy. It is lightweight, easy to operate, and ideal for patients who need reliable oxygen support while staying active and mobile.

Buy Now



HEMOCARE EQUIPMENT

OPEP Device

An OPEP device to help loosen and clear mucus from the airways, improving breathing, reducing congestion, and supporting better lung function.

Blood Pressure Monitors

A blood pressure monitor is used to regularly track a patients' blood pressure at home, helping detect changes early, manage chronic conditions, and support effective treatment.

Nebuliser

A nebuliser delivers medication directly to the lungs in a mist form, making it easier to treat respiratory conditions like asthma, COPD, or bronchitis and improve breathing.

Suction Unit

A home suction unit safely removes mucus and secretions from the airway, helping maintain clear breathing and prevent respiratory complications at home.

Pulse Oximeter

A pulse oximeter is used to quickly and non-invasively monitor their blood oxygen levels and heart rate, helping detect low oxygen and manage respiratory or cardiovascular conditions effectively.



OPEP Device

The Trudell Aerobika is a handheld oscillating PEP device designed to help clear mucus from the airways.

It uses gentle vibrations and positive pressure during exhalation to loosen secretions, making breathing easier and supporting effective airway clearance for patients with respiratory conditions.



Trudell Aerobika OPEP Device



Drug Free



Adjustable Settings



No Side Effects



Reusable

This clinically proven airway clearance device uses oscillating PEP therapy to help loosen and mobilize mucus for easier expectoration, features adjustable resistance levels to allow personalized therapy tailored to different patient needs and lung capacities, is drug-free and easy to use with a simple handheld design that supports effective daily therapy at home or in clinical settings, works in any position—sitting, standing, or lying down—to improve patient compliance, and is built to be durable and easy to clean, with simple disassembly for proper hygiene and repeated use.

Buy Now




How the *Aerobika*® device works.

The *Aerobika*® device has been clinically proven to improve ventilation, lung function and quality of life as well as reduce COPD flare-ups.¹




BEFORE USE: Airways inflamed, collapsed and plugged with mucus



AFTER USE: As you exhale through the device, the airways are opened



As the airways open, vibrations help thin and loosen mucus which can



This may help your inhaled medication work better²

Spacer Devices

The Trudell AeroChamber improves medication delivery to the lungs, reduces the need for hand-breath coordination, and features an anti-static chamber for consistent dosing, a built-in flow indicator for correct inhalation, and compatibility with masks, making it easy to use, hygienic, and highly effective for patients of all ages.



Trudell Aerochamber

CHILD



The Trudell AeroChamber is a valved holding chamber (VHC) designed to work with metered-dose inhalers (MDIs) to make inhaled medication more effective and easier to use. When a patient actuates their inhaler into the AeroChamber, the device holds the medication in a chamber, allowing the patient to inhale it slowly and steadily. This reduces the amount of medicine that ends up in the mouth or throat, ensuring more reaches the lungs. Its one-way valve prevents exhaled air from pushing the medicine back, while the anti-static chamber stops medication from sticking to the walls, so each dose is consistent.

Many models also come with mask attachments, making it suitable for children or patients who cannot use a mouthpiece properly. Overall, the AeroChamber simplifies inhaler use, improves treatment effectiveness, and ensures patients get the full benefit of their respiratory medication.



Drug
Free



Adjustable
Settings



No Side
Effects



Reusable

Buy Now



Spacer Devices

The Trudell AeroChamber improves medication delivery to the lungs, reduces the need for hand-breath coordination, and features an anti-static chamber for consistent dosing, a built-in flow indicator for correct inhalation, and compatibility with masks, making it easy to use, hygienic, and highly effective for patients of all ages.



Trudell Aerochamber



The Trudell AeroChamber is a valved holding chamber (VHC) designed to work with metered-dose inhalers (MDIs) to make inhaled medication more effective and easier to use. When a patient actuates their inhaler into the AeroChamber, the device holds the medication in a chamber, allowing the patient to inhale it slowly and steadily. This reduces the amount of medicine that ends up in the mouth or throat, ensuring more reaches the lungs. Its one-way valve prevents exhaled air from pushing the medicine back, while the anti-static chamber stops medication from sticking to the walls, so each dose is consistent.

Many models also come with mask attachments, making it suitable for children or patients who cannot use a mouthpiece properly. Overall, the AeroChamber simplifies inhaler use, improves treatment effectiveness, and ensures patients get the full benefit of their respiratory medication.



Drug
Free



Adjustable
Settings



No Side
Effects



Reusable

Buy Now



Blood Pressure Monitors

Ecomed Medical's blood pressure monitors provide accurate, reliable readings anytime, anywhere. Compact, easy-to-use, and clinically validated, they help you keep track of your health—backed by our nationwide service network for support and peace of mind.



Yuwell Upper Arm Blood Pressure Monitor



One Touch
Operation



Memory
Storage



Cuff
Detection



Irregular
Detection



Audio
Feedback

The Yuwell Upper Arm Blood Pressure Monitor provides accurate and reliable measurements of blood pressure and pulse, featuring an easy-to-read LCD display, one-touch automatic operation, memory storage for multiple readings, irregular heartbeat detection, and a comfortable adjustable cuff, making it ideal for convenient daily home monitoring and effective health management.

Buy Now



1 Year Warranty

Nebulisers

Ecomed Medical's nebulisers provide safe, effective, and reliable respiratory therapy in the comfort of your home or wherever you need it.

Designed for ease of use, they are compact, lightweight, and quiet, making treatment simple for patients of all ages.



Yuwell Mesh Nebuliser



Dual
Power



Ultra-
Silent



Minimal
Waste



Multiple
Accessories

Featuring advanced mesh technology for fast and efficient medication delivery, this nebuliser is portable and lightweight for use anywhere, operates quietly for comfort even during sleep, minimises medication waste to maximise drug use, and offers flexible power options with battery or USB operation.

Buy Now



1 Year Warranty

Yuwell Piston Nebuliser



Deep Lung
Delivery



Ultra-
Silent



Minimal
Waste



Oil Free

This nebuliser delivers a fine mist that reaches deep into the lungs for effective treatment, features a reliable, oil-free compressor for consistent airflow, is compact, lightweight, and easy to use, operates quietly for comfort across all ages, and minimises medication waste to maximise drug use.

Buy Now



1 Year Warranty

AeroEclipses

Ecomed Medical's AeroEclipse nebulisers provide efficient, reliable respiratory therapy with breath-actuated technology designed to help optimise medication delivery while reducing waste.



Trudell AeroEclipses



Trudell AeroEclipse is a breath-actuated nebuliser designed to deliver aerosol medication more efficiently during respiratory therapy.

Unlike traditional nebulisers that continuously release medication, the AeroEclipse only produces aerosol when the patient inhales, helping to reduce medication waste and improve delivery to the lungs.

It is lightweight, portable, and easy to use, making it suitable for both hospital and homecare environments.

The AeroEclipse is commonly used in the treatment of respiratory conditions such as asthma, COPD, and bronchiectasis, helping patients receive effective and reliable nebulised therapy.



Reduce
Waste



Portable



Direct
Delivery



Reusable

Buy Now



Suction Units

Ecomed Medical's suction units deliver safe, efficient, and reliable airway clearance therapy wherever it's needed —whether at home, in a healthcare setting, or on the go.

Designed with powerful yet gentle suction, they help remove mucus and secretions effectively, promoting easier breathing and reducing the risk of respiratory complications.



Yuwell Battery Operated Suction



Strong Suction



Dual Power



Safety Features



1000ml Jar



Ultra-Silent

Featuring dual power operation on battery or AC/DC, this suction unit delivers powerful suction to efficiently remove mucus and secretions, is portable and durable with a lightweight design and built-in handle, incorporates hygienic overflow protection and bacterial filters, and operates quietly with simple, easy-to-use controls for comfortable use.

Buy Now



2 Year Warranty

Yuwell Electric Suction



Strong Suction



Oil Free



Safety Features



1000ml Jar



Ultra-Silent

This suction unit delivers powerful, adjustable suction for effective mucus removal, features a maintenance-free, oil-free pump for reliable and easy care, includes safety features such as overflow protection and bacterial filters, offers user-friendly front panel controls, and combines quiet operation with a durable, long-lasting design.

Buy Now



2 Year Warranty

Pulse Oximeter

Ecomed Medical's pulse oximeters provide fast, accurate, and reliable monitoring of blood oxygen levels and heart rate —anywhere, at home or on the go.

Compact, lightweight, and easy to use, they allow patients and caregivers to track vital signs with confidence, helping detect changes early and supporting better overall health management.



Biolight M70 Finger Tip Pulse Oximeter



-  One Step Operation
-  Auto Power Off
-  Battery Powered
-  Accurate

This pulse oximeter provides fast and accurate SpO₂ and pulse readings, features a clear OLED display that's easy to read even in low light, includes auto power-off to save battery life, is portable and lightweight for home or travel use, and is built with a durable, shock-proof design for long-lasting reliability.

Buy Now



1 Year Warranty

Biolight M800 Handheld Pulse Oximeter



-  Baby
-  Dual Monitoring
-  Visual & Audio Alerts
-  Memory Storage

This versatile monitoring device measures SpO₂, pulse, ECG, and respiratory rate, features a large, easy-to-read colour OLED display, supports both spot checks and continuous monitoring, offers long battery life of up to 16 hours, and is portable and lightweight for convenient use at home, in clinics, or on the go.

View Brochure



1 Year Warranty

GAUTENG

Head Office - Krugersdorp

25 Clew Street, Monument Ext.1, 1745
Tel: 011 955 5710 / On-Call: 082 729 0749
Office Hours: Mon- Fri 08:00 - 16:30

Master Warehouse

Unit 5, Block B, Factoria Industrial Estate,
25 Louis Friedman Street, Factoria, 1739
Mon - Fri 08:00 - 16:30, Sat - Sun 08:00 - 12:00

Homecare Technical

Unit 9, Block C, Factoria Industrial Estate,
25 Louis Friedman Street, Factoria, 1739
Mon- Fri 08:00 - 16:30

East Rand

Unit 23, Oakhurst Office Park, 5 Grippen Rd, 1459
Tel: 011 027 4259
Mon - Fri 08:00 - 16:30, Sat 08:00 - 13:00

Pretoria Moot

271 15th Avenue, Rietfontein, 0084
Tel: 012 993 0493
Mon - Fri 08:00 - 16:30, Sat 08:00 - 13:00

Pretoria East

Unit 35, Scientia Techno Park,
2 Meiring Naudé Rd, Lynnwood, 0157
Tel: 012 885 0890
Mon- Fri 08:00 - 16:30

Vaal

137 General Hertzog Rd, Three Rivers,
Vereeniging, 1935
Tel: 016 013 0090
Mon- Fri 08:00 - 16:30

Centurion

13 Esdoring Nook, Highveld Technopark,
Centurion, 0157
Tel: 086 111 4028
Mon- Fri 08:00 - 16:30

FREE STATE

Bloemfontein

9 Nig Mientjie Street, Pellissier, 9301
Tel: 051 422 4118
Mon- Fri 08:00 - 16:30

Welkom

144 Constantia Street, Dagbreek, 9459
Tel: 057 030 0130
Mon- Fri 08:00 - 16:30

Bethlehem

43 De Leeuw St., Hospitaalheuwel, 9701
Tel: 082 593 8068
Mon- Fri 08:00 - 13:00

MPUMALANGA

Nelspruit

Unit 10 Central Park, Suikerriet Street, 1200
Tel: 010 446 9832
Mon- Fri 08:00 - 16:30

Witbank

1 Maria Street, Fransville, 1035
Tel: 013 590 0559
Mon- Fri 08:00 - 16:30

KWA-ZULU NATAL

Durban

450A Che Guevara, Glenwood, 4083
Tel: 031 202 4214
Mon- Fri 08:00 - 16:30

Dundee

84 Victoria Street, 3000
Tel: 082 491 3515
Mon- Fri 08:00 - 16:30

Pietermaritzburg

296 Bulwer St, CBD, 3201
Tel: 067 769 9057
Mon- Fri 08:00 - 13:00

WESTERN CAPE

Milnerton

Unit 36 Frazzitta Business Park, Freedom Way,
Milnerton, Cape Town, 7441
Tel: 021 020 1750
Mon- Fri 08:00 - 16:30

Muizenburg

4 Tokai Close, Marina de Gama, 7945
Tel: 021 788 3163
Mon- Fri 08:00 - 16:30

George

22 Waboom St, Glen Barrie, George, 6529
Tel: 044 124 0160
Mon- Fri 08:00 - 16:30

NORTH WEST

Klerksdorp

Unit 7, Stargate Office Park,
20 Buffeldoorn Avenue, Wilkoppies, 2571
Tel: 010 446 9833
Mon- Fri 08:00 - 16:30

Rustenburg

Unit 4, 250 Beyers Naude Dr,
Rustenburg, 2999
Tel: 014 112 0043
Mon- Fri 08:00 - 16:30

EASTERN CAPE

Port Elizabeth

29 Starling Crescent, Greenshields Park,
Gqeberha, 6070
Tel: 086 111 1718
Mon- Fri 08:00 - 16:30

East London

75 Devereux Avenue, Vincent, 5241
Tel: 010 446 9834
Mon- Fri 08:00 - 16:30

LIMPOPO

Polokwane

56 Bendor Drive, Bendor, 0699
Tel: 010 446 9831
Mon- Fri 08:00 - 16:30


NORTHERN CAPE

Kimberley

14A Schmidtdrift Rd, Verwoerdpark, 8301
Tel: 072 141 2335
Mon- Fri 08:00 - 16:30

Upington

25 Mark Street, Upington, 8801
Tel: 068 465 6688
Mon- Fri 08:00 - 16:30



Ecomed

Medical (Pty) Ltd



Through our 24 branches nation wide, we ensure patients have continuous, uninterrupted therapy.



+27 11 955 5710

We offer a 24-hour helpline, ensuring a patient's needs are supported at any time, day or night.

Visit our website



www.ecomed.co.za