

INFUSION AND BLOOD WARMER

MODEL: BFW-1000



BFW-1000

PRODUCT DESCRIPTION

- To warm up blood, blood products, pharmaceutical solutions, nutrition solutions or flushing liquid that are under refrigeration or at room temperature to the ideal temperature before infusion into the body to prevent hypothermia and its related complications.
- Blood and fluid transfusions for preoperative, intraoperative and postoperative patients; intravenous nutrient infusions; infusions for children or newborns; blood and fluid transfusions in cold environments.

PRODUCT FEATURES

Temperature control technology: the use of intelligent microcomputer control technology, triple temperature control protection

Dry heating structure, which heats up to the set temperature in 2 minutes

Alarm function: the infusion is finished and no utility power, the acousto-optic alarm to remind

Setting range: 28°C- 41°C, $\pm 0.5^{\circ}\text{C}$.

Heating capacity: 0-900 ml/H

BFW-1000: using high brightness LED digital tube to clearly display the liquid heating temperature

Fault acousto-optic alarm, built-in battery, infrared detection function can still work continuously for 3-6 hours after external power failure

TECHNICAL PARAMETERS

Executive Standard: YZB/Yue 1576-2012

Safety classification: Class I BF type

Electricity: AC 220V 50Hz Power: 70VA

Electric shock protection category: BF

Waterproof and dustproof rating: IP64

Fuse specification parameters: 2 X 1AL 250V

Battery parameters: 7.2V/500mAh NiMH rechargeable battery

Temperature setting range: 28°C-41°C (in 0.1°C increments)/82°F-105.8°F (in 0.2°F increments)

Measurable temperature range: -5~70°C

Working mode: continuous operation

Working Environment:

Temperature: +5°C~+25°C (41°F~104°F) Humidity: $\leq 80\%$ Atmospheric pressure: 86kPa~106kPa

Transport and storage environment

Temperature: -20°C~+55°C (14°F~104°F) Humidity: $\leq 80\%$ Atmospheric pressure: 86kPa~106kPa

No corrosive gas, good ventilation.



BFW-1000 LED Display

Shenzhen Bestman Instrument Co.,Ltd.

Manufacturer: Shenzhen Bestman Instrument Co.,Ltd.

Address: 8/F, Yifang Building, No. 315 Shuangming Avenue, Dongzhou Community, Guangming Street, Guangming District, Shenzhen, China.

Tel: 0755-26416184/26411423

Fax: 0755-26416184/26411423 Postcode: 518107

Web: www.szbestman.cn E-mail: bestman@szbestman.com



Wechat Public Account