

INFUSION AND BLOOD WARMER

MODEL: BFW-1020



BFW-1020

PRODUCT DESCRIPTION

- To warm up blood, blood products, pharmaceutical solutions, nutrition solutions or flushing liquid that are under refrigeration or at room temperature to the ideal temperature before infusion into the body to prevent hypothermia and its related complications.
- Blood and fluid transfusions for preoperative, intraoperative and postoperative patients; intravenous nutrient infusions; infusions for children or newborns; blood and fluid transfusions in cold environments.

PRODUCT FEATURES

- Temperature setting range of 28-41 degrees Celsius (in 0.1 increments)
- Fast: 2 minutes to set temperature, 24 hours continuous operation
- Equipped with high-flow transfusion extension special tube, which can increase the heating capacity in large quantity and speed.
- High efficiency: Intelligent microcomputer can automatically thermostat control, heating capacity up to 9 L/h, can simultaneously heating multiple isothermal liquid Unique drip rate/room temperature/function design.
- 3.2-inch color LCD display, real-time temperature, set temperature, heating time, alarm indication, room temperature, drip rate;
- Indication of with mains electricity or not, of with liquid or not and of with communication or not.

TECHNICAL PARAMETERS

Executive Standard: YZB/Yue 1576-2012

Safety classification: Class I BF type

Electricity: AC 220V 50Hz Power: 300VA

Electrocution protection category: BF Waterproof and dustproof rating: IP64

Fuse specification parameter: 2 X 1AL 250V

Battery parameters: 8.4V/220mAh NiMH rechargeable battery

Temperature setting range: 28°C-41°C (in 0.1°C increments)/82°F-105.8°F (in 0.2°F increments)

Working mode: Continuous operation

Operating environment:

Temperature: +5°C~+25°C (41°F~77°F) Humidity: ≤80% Atmospheric pressure: 86kPa~106kPa

Transportation and storage environment

Temperature: -20°C~+55°C (4°F~104°F) Humidity: ≤80% Atmospheric pressure: 86kPa~106kPa

No corrosive gas, good ventilation



3.2" color LCD Display

Shenzhen Bestman Instrument Co.,Ltd.

Manufacturer: Shenzhen Bestman Instrument Co.

Address: 8/F, Yifang Building, No. 315 Shuangming Avenue, Dongzhou Community, Guangming Street, Guangming District, Shenzhen, China.

Tel: 0755-26416184/26411423

Fax: 0755-26416184/26411423 Postcode: 518107

Web: www.szbestman.cn E-mail: bestman@szbestman.com



Wechat Public Account