



JX820D-1 Portable Aspirator(AC/DC)

Shanghai Medical Device Registration No. 20192140499

This model is the latest portable AC/DC suction unit. It is specially used for secretion, blood, vomitus and mucus removal of the patient who have undergone tracheostomy or had the trachea blocked, and for emergency treatment to a respiratory tract patient in hospital, first-aid teams, or social medical services

Product overview

1. Long life piston cylinder which requires no lubrication or daily maintenance produces no pollution to the environment.
2. Provided with three kinds of electrical sources, AC power, external DC power, and internal battery. They may be attached to the cigar lighter (DC12V) directly when used on the vehicles such as ambulance.
3. The internal battery is charged with constant voltage and limited current while the unit is working or being attached to AC100V~240V 50/60Hz or external DC 12V power supply. The capacity of the battery is respectively indicated by green, yellow and orange LEDs.
4. Overflow protection device can safely and efficiently prevent the liquid from entering into pump. Bottles made of plastic material are transparent, durable and convenient for cleaning.
5. Body made of plastic material is handy and beautiful. The wall-mounting frame makes it widely used in the room, on vehicle and conveniently hung on the side of the wheelchair.

Specification parameter

1. Max Vacuum: $\geq 0.08\text{MPa}$ (600mmHg)
2. 0.01Mpa(75mmHg~Max vacuum)
3. Flow rate: $\geq 20\text{ L/min}$
4. Noise: $\leq 65\text{dB (A)}$
5. Jar capacity: 1000mL (PC)
6. Operating voltage: DC 12V ; AC 100V~240V 50/60Hz
7. Input: 110VA
8. Lithium battery, Battery working time (duration on full charge): $\geq 60\text{min}$
9. Dimension: 350mm x 140mm x 340mm
10. Package Dimension: 405mmx200mmx360mm
11. G.W./N.W.: 5/4.5kg



(Optional)

SHANGHAI BAOJIA MEDICAL APPARATUS CO.,LTD.

Add.: 139 Yuantai Road, Baoshan District, Shanghai, China

Zip code: 200444

Tel: +86 21-66740008

E-mail: sales02@smaf.com.cn

Website: www.smaf.com.cn

