



KL-LED MSTZ4 AC/DC
Shadowless Operating Lamp

KL-LED•MSTZ4 AC/DC SHADOWLESS OPERATING LAMP

KL-LED.MSTZ4 AC/DC LED shadowless operating lamp is suitable for various operation needs. It is a modern ideal surgical lighting equipment operating room and purify the operating room.

Features ●

- Adopt imported LED cold light source, white the ultra-thin optical lens
- Adopt digital means for stepless regulation of LED brightness.
- Adjust the diameter of light spot and color temperature
- With Endoscopy model
- Removable handle sheath, easy to clean and disinfect
- Nano coating of heat sink has created an excellent heat dissipation effect.
- Imported balance arm, which is flexible up and down

Specifications ●

Mode	KL-LED.MSTZ4 AC/DC
Illumination at 1m(lux)	40000~160000
Light field diameter at 1m (cm)	60-100cm
D10(mm)	180±40
Diameter of spot D50(mm)	110±35
Color rendering index(Ra)	100≥Ra≥95
Optional color temperature(K)	6700≥Tc≥3000K
Endo model(lux)	10000
Number of LED bulbs	45
Average bulb life (hours)	50,000
Heat-to-light ratio(mW/m ² lx)	<6
Working distance(cm)	70-140
Max. power consumption(VA)	90



● KL- LED•MSTZ4 (AC/DC)



让爱照亮未来
Future Is Well Lit By Love

上海科凌国际贸易有限公司

Shanghai Keling International Trading Co., Ltd.

Company Address: No. 1725 Huangxing Road, Yangpu District, Shanghai, China

Factory Address:

Factory 1: Shanian Industrial Zone, Gaoyou City, Yangzhou City, Jiangsu Province, China

Factory 2: No. 18, Xiaoguanqiao Road, Jiangyang Industrial Zone, Yangzhou City, Jiangsu Province, China

Email: webmaster@shkeling.com

Mob: +86 15000649020

Web: en.jskeling.com



Keling Website



Keling Facebook