





KL-LED.SZ4/STZ5
Shadowless Operating Lamp

KL-LED.SZ4/STZ5 Shadowless Operating Lamp

Instructions

◆ KL-LED.SZ4/ STZ5 LED shadowless operating lamp is suitable for various operation needs.
It is a modern ideal surgical lighting equipment operating room and purify the operating room.



Mode	KL-LED.SZ4/STZ5
Illumination at 1m(lux)	40000~120000 / 40000~160000
Illumination at 1m on Endo mode(lux)	--/ 10000
Light field diameter at 1m (cm)	10-28
Light filed diameter D10mm	180±40 / 180±40
D50(mm)	110±35 / 110±35
Depth of L1+L2 illumination(cm)	150
Color rendering index(Ra)	100≥Ra≥85
Optional color temperature(K)	4500±500K 5000≥Tc>3000K
Number of LED bulbs	32/45
Average bulb life (hours)	50,000
Heat-to-light ratio(mW/m ² lx)	<6
Working distance(cm)	70-140
Power supply	100-240 V, 50/60HZ
Max. power consumption(VA)	80W / 90W
Ceiling height for installation(cm)	270-320



让爱照亮未来
Future Is Well Lit By Love

上海科凌国际贸易有限公司

Shanghai Keling International Trading Company

Company Address: No. 1725 Huangxing Road, Yangpu District, Shanghai, China

Factory Address:

Factory 1: Shanian Industrial Zone, Gaoyou City, Yangzhou City, Jiangsu Province, China

Factory 2: No. 18, Xiaoguanqiao Road, Jiangyang Industrial Zone, Yangzhou City, Jiangsu Province, China

Email: webmaster@shkeling.com

Mob: +86 15000649020

Web: en.jskeling.com



Keling Website



Keling Facebook