

O2FLO PRO™

High Flow Heated Respiratory Humidifier

Model No.: VUN-003-A

The total solution 80LPM therapy including a complete range of consumables

- **Electronic Auto-blending** for
 - ✓ Easy setup
 - ✓ Accurate FiO₂ delivery
- Real time **SpO2 detection (optional)**
 - ✓ **ROX index** as indication of potential treatment effectiveness
- Proprietary **infrared water level detection** allowing time to refill, and PREVENTS
 - ✗ dry-out
 - ✗ overheating
- 4800mAh Battery for **Transport Mode**
- Continuous respiratory support during short distance transport
- **Low Flow Mode** for pre-set treatment interval settings
- Suitable for THRIVE, post-op, patient weaning, COPD, and patients suffering from any other kinds of ARDS
- Interactive touchscreen
- **HEPA-type filter** for clean breathing
- **USB** port for extended function
- **SD card** to record alarm and usage memory

Safer · Smarter

The All-New Nasal Cannulas for Better Comfort & Performance

- Heated Breathing Circuits and Nasal Cannula that delivers up to 80LPM
- Optimal humidification and high flow results with inspired™ package





ADULT

NASAL CANNULA	SIZE	Maximum flow rates(LPM)
● VANC-01	S	Up to 75
● VANC-02	M	Up to 80
● VANC-03	L	Up to 80

- Increased maximum flowrate up to 80LPM
- Soft cheek pads for extra comfort
- 270°- swivel allowing customised fitting for all cheek shapes



PAEDIATRIC / INFANT

NASAL CANNULA	SIZE	Maximum flow rates(LPM)
● VINC-05	XS	Up to 8
● VINC-06	S	Up to 9
● VINC-07	M	Up to 10
● VINC-08	L	Up to 20
● VINC-09	XL	Up to 25

- Increased maximum flowrate up to 25LPM
- Color coded for immediate sizing visualization
- Soft and large cheekpads for reducing pressure on face

Accessories



Adult Nasal Cannula
VANC-01 (S)
VANC-02 (M)
VANC-03 (L)



Adult Nasal Cannula (Premium)
● VANC-01 (S)
● VANC-02 (M)
● VANC-03 (L)



Infant Nasal Cannula
○ VINC-01 (S) ● VINC-02 (M)
● VINC-03 (L) ● VINC-04 (XL)



Infant Nasal Cannula (Premium)
● VINC-05 (XS) ● VINC-06 (S)
● VINC-07 (M) ● VINC-08 (L)
● VINC-09 (XL)



Tracheostomy Interface
VTI-001



Integrated Heated Breathing Circuit
Audit: 510-049
Paediatric: 510-050



Embedded Heated Breathing Circuit
Audit: 510-049S
Paediatric: 510-050S



Double Wall Heated Breathing Circuit
Audit: 510-049D
Paediatric: 510-050D



**Flow Meter BS
Flow Meter GS
Flow Meter DISS**



Oxygen Tube



Humidification Chamber
VHC60



Double Wall Humidification Chamber
VHC65



Humidification Chamber
VHC80
(Can be Used with Nebulizer)



HEPA-type Filter



Connection Cable
VUN-300



Trolley with Arm
MT001

*Product availability may vary amongst territories.

Contact us for more information

email: info@inspired-medical.com

www.vincentmedical.com

www.inspired-medical.com

REF# IMDM004 Rev.A 2024-09



Vincent Medical



inspired™ Medical