



NEBTIME UN600A



Product Page

Evaporation of liquids and medicine in solution form allows it to be used in the treatment of patients.



Digital Display

Thanks to the LCD screen, 4 main treatment settings such as aerosol, air flow, heater and timer can be easily set.



Feeding Bottle

Thanks to its 1 liter capacity jar, it extends the duration of the treatment and provides a more effective use.



Autoclavable Hose

It provides a more sterile use, thanks to the autoclavable hose that withstands 121 degrees Celsius.



Hospital Use

Designed with hospital staff in mind, Nebtime UN600A offers a comfortable use in the hospital.



Smooth Design

Elmaslar Medical Devices has won numerous design awards for its ergonomic design approach, which is evident in all of its products.

Technical Specifications

Voltage	230 VAC $\pm$ 10%, 50 Hz/60Hz
Power	40 W - 70 W
Crystal Operating Frequency	1,68 MHz $\pm$ 5%
Particle Size	0,5 – 6 μ m
Average Nebulization Rate	0-5 ml/dk $\pm$ 10%
Aerosol Tank	650 ml
Weight	2.75 kg
Gross Weight	11 kg
Heater Max.	38°C (with heater)
Operating Temperature Range	10°C ~ 40°C
Storage Temperature Range	10°C ~ 40°C
Power Consumption	70W
Operation Type	Continuous
Device Dimensions	27x23x21 cm
Box Dimensions	64x62x23 cm
Stock Code	UN01B000
Warranty Period	2 Years
Noise Level	< 40 dB
Medicine Cup Capacity	70 ml
Maximum Medicine Capacity	30 ml