



BREATHCARE ECO II

ET5

AUTO CPAP



Working Mode:	CPAP, Auto CPAP
Pressure Range (cm H ₂ O):	4 to 20
Humidifier capacity:	300ml
Noise level:	28dB(A)
Net weight:	1.4kg
Dimension (cm):	27x13x14
Certification:	CE
Warranty:	Device (3 years) Water tank and air Tubing (90days)

Auto-All System	Smart Start/Stop Function
Enhanced Adaptive Algorithm	Leak Compensation
Wireless data system	Remote Monitoring

CONFIGURATION



Water Tank



Power Adapter



Air Tubing



SD Card



Air Filter (2 pcs)



Portable Bag



Free Platform

OPTIONAL

- 1: Transmission module (WiFi function/ 3 yrs SIM)
- 2: Heated tubing
- 3: BreathWear Mask



New Edition, Awesome Value

Next Industrial Design

Built on a next-gen ECO design logic, this system ensures accurate pressure delivery while maintaining a precisely silent operation. It harmonizes advanced performance stability with intuitive interaction, providing a consistently reliable and comfortable therapy experience that promotes daily adherence.



Enhanced Adaptive Algorithm

Utilizes advanced signal processing to maintain precise pressure control across dynamic respiratory events.



Intuitive Auto-All Therapy

A seamless fusion of automated intelligence and user-centered design. The system manages all therapy adjustments in the background, delivering a "set-and-forget" experience that is both simple to use and easier to stick with, ultimately improving treatment adherence.

