

Examination light, LED, Mobile, with Battery

KL01L.IILED (AC/DC) Suitable for ENT, urology, gynecology examination and outpatient surgery half auxiliary lighting.



- ◆ An imported LED cold light source is adopted as operation lighting. As a genuine cold light source.
- ◆ The LED shadowless lamp is pure-DC power supplied, ensuring no flickering or dimming occurs.
- ◆ The service life of LED shadowless lamp is long (50,000 h).
- ◆ LED has high luminous efficiency and it is impact-resistant and strong, without mercury contamination. The light it sends contains no infra-red or ultraviolet radiation pollution.
- ◆ The unique design of the spring arm. Allows for easy, single-handed manoeuvring.
- ◆ A weighted base allows for manoeuvrability without compromising on balance.

Technical Specifications

Effective lighting over the bed with a minimum length of 1m

The light may be separately switched on via the column or the luminaire.

Colour rendering index: > 90

Working distance of 50 cm = 15 cm

Working distance of 100 cm = 25 cm

Operates from 230V, 50 Hz mains power

Mobile stand: lockable castors + hook for mains cable

Durable castors

Smooth surface allows for easy cleaning

Colour temperature 4000k

Lamp life:>50000 Continuous hours,

Base diameter 40cm, Battery life 4 hours

60000 lux

