

Examination light, LED, Wall Mount

KL01L.III LED Wall Mount suitable for ENT, urology, gynecology examination and outpatient surgery half auxiliary lighting.

- Features**
- ◆ An imported LED cold light source is adopted as operation lighting. As a genuine cold light source.
 - ◆ The LED shadowless lamp is pure-DC power supplied, ensuring no flickering or dimming occurs.
 - ◆ The service life of LED shadowless lamp is long (50,000 h).
 - ◆ LED has high luminous efficiency and it is impact-resistant and strong, without mercury contamination. The light it sends contains no infra-red or ultraviolet radiation pollution.
 - ◆ The unique design of the spring arm. Allows for easy, single-handed manoeuvring.
 - ◆ Unique lamp design allows for efficient, out of the way storage.

Technical Specifications

Effective lighting over the bed with a minimum length of 1m. Colour rendering index: > 90

Working distance of 50 cm = 15 cm

Working distance of 100 cm = 25 cm

Operates from 230V, 50 Hz mains power, built in power supply

Smooth surface allows for easy cleaning

Colour temperature 4000k

Lamp life : >50000 Continuous hours

60000 lux



Sterilisable positioning handle